




Entrada a Higuito • San Isidro del Guarco, 30802 • Cartago, Costa Rica
800-278-9947 • mail@glidewell dental.cr • cr.glidewell dental.com

Dr. Name _____ Phone # _____
Patient Name _____ Acct. # _____
Address _____
Email _____ Deliver by 5 p.m. on _____
Enclosed with Case: Impressions Models Bite Photos Other: _____ Test: Yes No

SHADE INSTRUCTIONS



Tooth No. _____
Final Shade _____
Stump Shade _____
(Required for BruxZir Esthetic and IPS e.max)

OCCLUSAL STAINING

None Light* Medium Dark

IF NO OCCLUSAL CLEARANCE

Call doctor Spot opposing
 Metal occlusion Metal island
 Make this a permanent note

PONTIC DESIGN



*

ALL-CERAMIC RESTORATIONS

BruxZir Full-Strength* (> 1,000 MPa)
 BruxZir Esthetic (≥ 870 MPa)
(Providing the stump shade is recommended for restorations less than 1.5 mm thick)
 IPS e.max
 Crown Bridge
 Inlay/Onlay Veneer

PORCELAIN FUSED TO METAL

Non-Precious Metal* Maryland bridge
 White Semi-Precious Metal

FULL CAST RESTORATION

Non-Precious Semi-Precious

INCLUSIVE CUSTOM ABUTMENTS

Titanium*
 Zirconia w/ Ti-Base
 Semi-Precious
 Prepare existing abutment

Implant System _____


Implant Diameter _____ mm

TCS FLEXIBLE REMOVABLES

iFlex
 Bilateral Unilateral
 Totally Natural
 Karadent
 All Pink Transparent/Pink
 Unbreakable
 Bilateral Unilateral
 Acrylic Shade _____
 Kenson Teeth
(Included at no extra charge) _____
 Premium brand teeth (Extra charge applies)
Shade _____ Brand _____ Mould _____
 Bite rim Wax setup try-in Finish

Rx SPECIFIC INSTRUCTIONS

NOTE: Please send a study model on all work involving anterior teeth.



SCREW-RETAINED RESTORATIONS

BruxZir Full-Strength* (w/Ti/Base)
 IPS e.max Press (w/Ti/Base)
 PFM Semi-Precious
 Temporary Screw-Retained

SELECT RESTORATION TYPE

Crown Splinted Crowns Bridge

TEMPORARIES

BioTemps Provisionals (PMMA)
 BioTemps with metal
 BioTemps Provisionals Screw Retained Crown
 Smile Transitions
 Diagnostic Wax-Up

SLEEP APPLIANCES

EMA Silent Nite

BITE SPLINTS

Upper* Lower
 Comfort H/S Bite Splint (hard/soft)
 Clear* Blue Pink
 Comfort Bite Splint (hard)
 Clear* Turquoise
 Clear Soft Nightguard
 Semi-Hard Nightguard
 Bleaching Tray
 Ortho Retainer

PLAYSAFE MOUTHGUARDS

Junior Light Light Pro
 Light Medium Heavy Heavy Pro
 Helmet strap _____

Specify color(s) and name on RX.

Upper Lower

Signature _____

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

License _____

Date _____

TERMS AND WARRANTY INFORMATION



Made with materials from the United States

We honor VISA, MASTERCARD, AMEX and DISCOVER.

TERMS: All accounts are payable within 15 days of statement date. Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance. Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Glidewell Costa Rica ("the lab") is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit cr.glidewelldental.com/sobre-nosotros/politicas-garantias.

PRIVACY POLICY: The data that is supplied will be treated in compliance with Law 8968 of Costa Rica. The names and patient data will be treated as a reference code and not as proper names or personal data. Details available to you on our website: cr.glidewelldental.com

SHIPPING: There is no charge for one inbound and one outbound shipment per case. Additional shipments for die trims, bisque bake, coping or framework try-ins, reshades to a new shade, and oversized articulators will incur additional fees.



• **BruxZir Restorations**



• **IPS e.max®**
• **PFM Restorations Obsidian®**




• **BioTemps® Provisionals (PMMA)**
• **Bite Splints**
• **Mouthguards and Nightguards**
• **tcs® Flexible Removables**



• **Custom Abutments**

All rush cases must be prescheduled by calling 800-278-9947 before the case is shipped. Time of pickup and delivery may affect turnaround time.

PREPARATION GUIDELINES



BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

FLAT-RATE PRICE ON CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS IS AVAILABLE FOR THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™
Certain®

CAMLOG®
SCREW-LINE

DENTSPLY Implants
ANKYLOS® C/X
ASTRA TECH Implant System®
ASTRA TECH Implant System® EV

Glidewell Direct
Hahn™ Tapered Implant System
Inclusive® Tapered Implant System

HIOSSEN®
HG System

MegaGen
AnyRidge® Implant System

Nobel Biocare
Brånemark System® RP
NobelActive®
NobelReplace®

Straumann®
Bone Level
Tissue Level

Zimmer Dental
Screw-Vent®

Restorations for the implant systems above are fabricated using components manufactured by Prismatic Dentalcraft, Inc. Prices may vary for other implant systems due to the need for original equipment manufacturer (OEM) components. Inclusive is a registered trademark of Prismatic Dentalcraft, Inc. Hahn Tapered Implant is a trademark of Prismatic Dentalcraft, Inc. All other trademarks are property of their respective owners.